



Bellingham+Stanley

Frequently Asked Questions

FAQ No: FAQ-GRP-056

Title: RoHS and WEEE European directives on lamps in category 5

Background:

This statement outlines the current position and policies of **Bellingham and Stanley Ltd** with regard to all products containing removable lamps sold in the European Community.

The following products manufactured and supplied by Bellingham and Stanley fall under the scope of the RoHS directive (Reduction of Hazardous Substances).

B & S Part No.	Description	No. of lamps
59-01	Free standing sodium lamp outfit 230v \pm 10% 50-60Hz	ONE
59-06	Free standing sodium lamp outfit 110v \pm 10% 50-60Hz	ONE
59-20	Free standing spectral lamp outfit 230v \pm 10% 50-60Hz	ONE
59-25	Free standing spectral lamp outfit 110v \pm 10% 50-60Hz	ONE
59-03	Model D base sodium lamp outfit 230v \pm 10% 50-60Hz	ONE
59-08	Model D base sodium lamp outfit 110v \pm 10% 50-60Hz	ONE

Note:

Removable lamps – the above table does not include potted LED's or LED's incorporated on a circuit board where the LED has been directly mounted. All our digital refractometers and polarimeters use potted LED's as a light source and directly mounted LED's for power indicators.